Fetal Cell Count™ kit

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Bibliography

1. NEN EN ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements.
2. Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.
3. Garratty, G., and Arndt, P.A. 1999. Applications of flow cytometry to red blood cell immunology. *Cytometry (Communications in Clinical Cytometry)* 38: 259-267.
4. Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.

5. Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.
6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf. Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin. Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.*48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetal maternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.
11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.
13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. *Vox Sang.* 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: *International Society of Blood Transfusion – 5th Regional Congress; Venice 2-5 July*, abstract POS 309, p239.
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretière, J., Guimbretière, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. *Rev.Fr.Transfus.Hémobiol.* 35: 239-254.
16. Brady, H.J.M., Edwards, M., Linch, D.C., 1990. Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. *British Journal of Haematology*, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athen, Griechenland, 20. -25. September 2005.
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. *Transfusion*, 47:7 , 1281 – 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. *Eur J Obstet Gynecol Reprod Biol.* 2006 Jul 24
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. *Am.J.Clin.Pathol.*136(4):631-6.

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Bibliografi

1. DIN EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.
2. Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.
3. Garratty, G., and Arndt, P.A. 1999. Applications of flow cytometry to red blood cell immunology. *Cytometry (Communications in Clinical Cytometry)* 38: 259-267.
4. Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.
5. Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.
6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf.Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin.Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.* 48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetalmaternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.
11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.

13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. Vox Sang. 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: International Society of Blood Transfusion – 5th Regional Congress; Venice 2-5 July, abstract POS 309, p239.
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretière, J., Guimbretière, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. Rev.Fr.Transfus.Hémobiol. 35: 239-254.
16. Brady, H.J.M., Edwards, M., Linch, D.C., 1990 Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. British Journal of Haematology, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athens, Greece, 20-25 September 2005
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. Transfusion, 47:7 , 1281 – 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. Eur J Obstet Gynecol Reprod Biol. 2006 Jul 24.
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. Am.J.Clin.Pathol.136(4):631-6.

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- DIN EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.
- Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.

3. Garratty, G., and Arndt, P.A. 1999. Applications of flow cytofluorometry to red blood cell immunology. *Cytometry (Communications in Clinical Cytometry)* 38: 259-267.
4. Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.
5. Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.
6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf.Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin.Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.*48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetomaternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.
11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.
13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. *Vox Sang.* 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: *International Society of Blood Transfusion - 5th Regional Congress; Venice 2-5 July, abstract POS 309, p239.*
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretière, J., Guimbretière, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. *Rev.Fr.Transfus.Hémobiol.* 35: 239-254.
16. Brady, H.J.M., Edwards, M., Lynch, D.C., 1990. Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. *British Journal of Haematology*, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athens, Greece, 20-25 September 2005.
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. *Transfusion*, 47:7, 1281 - 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. *Eur J Obstet Gynecol Reprod Biol.* 2006 Jul 24.
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. *Am.J.Clin.Pathol.*136(4):631-6

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©2021 - IQ Products bv. All rights reserved. No parts of these works may be reproduced in any form without permission in writing.

Current version + release date	Version 1 01-09-2021 (DD-MM-YYYY)
Previous version	N/A
Changes	1. Created this online version with translations of the international booklet which is printed and included with every kit. The online version includes Portuguese, Swedish and Spanish, whereas the international booklet includes English, French, German and Italian.
Justification	1. Environmental reasons.