# Fetal Cell Count™ kit

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## Bibliography

- EN ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements.
- Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.
- Garratty, G., and Arndt, P.A. 1999. Applications of flow cytometry to red blood cell immunology. *Cytometry (Communications in Clinical Cytometry)* 38: 259-267.
- Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.
- Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.

6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf. Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin. Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.* 48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetal maternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.
11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.
13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. *Vox Sang.* 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: *International Society of Blood Transfusion – 5<sup>th</sup> Regional Congress; Venice 2-5 July*, abstract POS 309, p239.
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretière, J., Guimbretière, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. *Rev.Fr.Transfus.Hémobiol.* 35: 239-254.
16. Brady, H.J.M., Edwards, M., Linch, D.C., 1990 Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. *British Journal of Haematology*, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athen, Griechenland, 20. -25. September 2005.
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. *Transfusion*, 47:7 , 1281 – 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. *Eur J Obstet Gynecol Reprod Biol.* 2006 Jul 24.
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. *Am.J.Clin.Pathol.* 136(4):631-6.

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## Bibliografi

1. EN ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements.
2. Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.
3. Garratty, G., and Arndt, P.A. 1999. Applications of flow cytometry to red blood cell immunology. *Cytometry (Communications in Clinical Cytometry)* 38: 259-267.
4. Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.
5. Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.
6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf.Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin.Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.* 48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetomaternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.

11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.
13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. *Vox Sang.* 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: *International Society of Blood Transfusion – 5<sup>th</sup> Regional Congress; Venice 2-5 July*, abstract POS 309, p239.
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretiére, J., Guimbretiére, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. *Rev.Fr.Transfus.Hémobiol.* 35: 239-254.
16. Brady, H.J.M., Edwards, M., Linch, D.C., 1990 Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. *British Journal of Haematology*, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athens, Greece, 20-25 September 2005
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. *Transfusion*, 47:7, 1281 – 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. *Eur J Obstet Gynecol Reprod Biol.* 2006 Jul 24.
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. *Am.J.Clin.Pathol.* 136(4):631-6.

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- EN ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements.
- Sebring, E.S., Polesky, H.F. 1990. Fetomaternal hemorrhage: incidence, risk factors, time of occurrence, and clinical effects. *Transfusion* 30: 344-357.
- Garratty, G., and Arndt, P.A. 1999. Applications of flow cytometry to red blood cell immunology.

Cytometry (Communications in Clinical Cytometry) 38: 259-267.

4. Nance, S.J., Nelson, J.M., Arndt, P.A., et al. 1989. Quantitation of Feto-maternal hemorrhage by flow cytometry, a simple and accurate method. *Am.J.Clin.Pathol.* 91: 288-292.
5. Hadley, A.G. 1998. A comparison of in vitro tests for predicting the severity of haemolytic disease of the fetus and newborn. *Vox Sang.* 74: 375-383.
6. Lee, D., Contreras, M., Robson, S.C., Rodeck, C.H., Whittle, M.J. 1999. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. *Transf.Med.* 9: 93-97.
7. Kleihauer, P., Braun, H., and Betke, K. 1957. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. *Klin.Wochenschr.* 35: 637-638.
8. Davis, B.H., Olsen, S., Bigelow, N.C., Chen, J.C. 1998. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. *Immunohematology* 38: 749-756.
9. Johnson, P.R., Tait, R.C., Austin, E.B., et al. 1995. Flow cytometry in diagnosis and management of large fetomaternal haemorrhage. *J.Clin.Pathol.* 48: 1005-1008.
10. Nelson, M., Popp, H., Horky, K., Forsyth, C., Gibson, J. 1994. Development of a flow cytometric test for the detection of D-positive fetal cells after fetomaternal hemorrhage and a survey of the prevalence in D-negative women. *Immunohematology* 10: 55-59.
11. Corsetti, J.P., Cox, C., Leary, J.F., Cox, M.T., Blumberg, N., Doherty, R.A. 1987. Comparison of quantitative acid-elution technique and flow cytometry for detecting fetomaternal hemorrhage. *Ann.Clin.Lab.Sci.* 17: 197-206.
12. Navenot, J.M., Merghoub, T., Ducrocq, R., Krishnamoorthy, R., Blanchard, D. 1998. A new method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle cell disease. *Cytometry* 32: 186-190.
13. Nelson, M., Zarkos, K., Popp, H., Gilson, J. 1998. A flow-cytometric equivalent of the Kleihauer test. *Vox Sang.* 75: 234-241.
14. Navenot, J.M., Blandin, A.M., Willem, C., Bernard, D., Muller, J.Y., Blanchard, D. 1995. In situ detection of fetal hemoglobin by flow cytometry: evaluation of a simple procedure for quantitating fetal erythrocytes in maternal peripheral blood. In: *International Society of Blood Transfusion – 5<sup>th</sup> Regional Congress; Venice 2-5 July*, abstract POS 309, p239.
15. Blanchard, D., Bernard, D., Loirat, M.J., Frioux, Y., Guimbretière, J., Guimbretière, L. 1992. Caractérisation d'anticorps monoclonaux murins dirigés contre les érythrocytes foetaux. *Rev.Fr.Transfus.Hémobiol.* 35: 239-254.
16. Brady, H.J.M., Edwards, M., Linch, D.C., 1990. Expression of the human carbonic anhydrase I gene is activated late in fetal erythroid development and regulated by stage-specific trans-acting factors. *British Journal of Haematology*, 76, 135-142.
17. Bernaud, J., Rigal, D., Porra, V., Follea, G., Blanchard, D. Fetal Cell Count™ - a commercial kit for quantification of fetal cells in maternal blood by flow cytometry. 5th Euroconference on Clinical Cell Analysis. Athens, Greece, 20-25 September 2005.
18. Porra, V., Bernaud, J., Gueret, P., Bricca, P., Rigal, D., Follea, G., and Blanchard, D. 2007 Identification and quantification of fetal red blood cells in maternal blood by a dual-color flow cytometric method: evaluation of the Fetal Cell Count kit. *Transfusion*, 47:7, 1281 – 1289.
19. Leers MP, Pelikan HM, Salemans TH, Giordano PC, Scharnhorst V. Discriminating fetomaternal hemorrhage from maternal HbF-containing erythrocytes by dual-parameter flow cytometry. *Eur J Obstet Gynecol Reprod Biol.* 2006 Jul 24.
20. de Wit H, Nabbe KC, Kooren JA, Adriaansen HJ, Roelandse-Koop EA, Schuitemaker JH, et al. 2011 Reference values of fetal erythrocytes in maternal blood during pregnancy established using flow cytometry. *Am.J.Clin.Pathol.* 136(4):631-6

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|                                |                                                                                                                                                                                                                                                                                                                                           |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Current version + release date | Version 4<br>01-11-2025 (DD-MM-YYYY)                                                                                                                                                                                                                                                                                                      |
| Previous version               | Version 3                                                                                                                                                                                                                                                                                                                                 |
| Changes                        | 1. The <b>CE</b> symbol is no longer shown in italics<br>2. hazard symbols of reagent B have been updated                                                                                                                                                                                                                                 |
| Justification                  | 1. In response to a non-conformity raised by the notified body.<br>2. In response to the CLP analysis performed for the IVDR Technical Documentation.                                                                                                                                                                                     |
| Current version + release date | Version 3<br>29-03-2023 (DD-MM-YYYY)                                                                                                                                                                                                                                                                                                      |
| Previous version               | Version 2                                                                                                                                                                                                                                                                                                                                 |
| Changes                        | 1. Added 'Regulatory status' and removed separate sentence about IVD registration.<br>2. Document title of the EN ISO 15223-1 has been updated.                                                                                                                                                                                           |
| Justification                  | 1. Better overview of IVD registrations.<br>2. Document title was adapted with current version.                                                                                                                                                                                                                                           |
| Current version + release date | Version 2<br>02-03-2022 (DD-MM-YYYY)                                                                                                                                                                                                                                                                                                      |
| Previous version               | Version 1                                                                                                                                                                                                                                                                                                                                 |
| Changes                        | 1. In the last sentence of paragraph 'Intended use' 2 hours replaced by 90 minutes (p.3, 8, 13),<br>2. In the paragraph 'Control sample preparation' the sentence 'The spike for the quality control of the kit is made with the appropriate volume of washed cells from 0.5 mL cord blood and adult blood.' has been added (p.4, 9, 14). |
| Justification                  | 1. The complete time for staining and analysis is 90 minutes, as already mentioned on the website and the product flyer.<br>2. The sentence has been added as advice for the users at the laboratory.                                                                                                                                     |
| Current version + release date | Version 1<br>01-09-2021 (DD-MM-YYYY)                                                                                                                                                                                                                                                                                                      |
| Previous version               | N/A                                                                                                                                                                                                                                                                                                                                       |
| Changes                        | Created this online version with translations of the international booklet which is printed and included with every kit. The online version includes Portuguese, Spanish and Finnish, whereas the international booklet includes English, French and German.                                                                              |
| Justification                  | Environmental reasons.                                                                                                                                                                                                                                                                                                                    |

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©2025 - IQ Products bv. Todos os direitos reservados Os elementos contidos neste documento não podem ser reproduzidos sem consentimento por escrito.