FMH Quikquant™

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Referências

- 1** CLSI, Fetal red cell detection; approved guideline. CLSI (formerly NCCLS) Document H52-A, 2001.
- 2** Davis, BH. Diagnostic advances in defining erythropoietic abnormalities and red cell diseases. Seminars in Hematology, 2001;38:148-59.
- 3** Davis, B.H. Diagnostic utility of red cell flow cytometric analysis. Clin Lab Med 2001;21(4):829-40.
- 4** Sebring ES, Polesky HF. Fetomaternal haemorrhage: incidence, risk factors, time of occurrence, and clinical effects. Transfusion 1990;30:344-57.
- 5** Giacola G.P. Severe fetomaternal hemorrhage: a review. Obstet & Gynecol Surv, 1997;52(6): p. 372-80.
- 6** Polesky, HF, Sebring ES, Evaluation of methods for detection and quantitation of fetal cells and their effect on RhIgG. American Journal of Clin Path, 1981;76:525-29.
- 7** Hartwell, EA. Use of Rh immune globulin: ASCP practice parameter. American Journal of Clin Path, 1998;110:281-92.
- 8** Lee D, Contreras M, Robson SC, et al. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. Transfusion Med 1999;9:93-7.
- 9** Kleihauer E, Braun H, Betke K. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. Klin Wochenschr 1957;35:637-38.
- 10** Duckett JR, Constantine G. The Kleihauer technique: an accurate method of quantifying fetomaternal haemorrhage? British Journal of Obstet & Gynaecol 1997;104: 845-6.
- 11** Emery CL, Morway LF , et al. The Kleihauer-Betke test. Clinical utility, indication, and correlation in patients with placental abruption and cocaine use. Arch Pathol Lab Med 1995;119(11):1032-7.
- 12** Davis BH, Olsen S, et al. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. Transfusion 1998;38(8):749-56.
- 13** Chen JC, Davis BH, et al. Multicenter clinical experience with flow cytometric method for fetomaternal hemorrhage detection. Cytometry 2002;50(6):285-90.
- 14** Lloyd-Evans P, Kumpel BM et al. Use of a directly conjugated monoclonal anti-D (BRAD-3) for quantification of fetomaternal hemorrhage by flow cytometry. Transfusion 1996;36(5):432-7.
- 15** Nance SJ, Nelson JM, et al. Quantitation of fetal-maternal hemorrhage by flow cytometry. A simple and accurate method. Am J Clin Pathol, 1989;91(3): 288-92.
- 16** Navenot JM, Merghoub T et al. New method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle-cell disease. Cytometry 1998;32(3):186-90.
- 17** Nelson M, Zarkos K et al. A flow-cytometric equivalent of the Kleihauer test. Vox Sanguinis 1998;75:234-41.
- 18** Navenot JM, Muller JY, et al. Expression of blood group i antigen and fetal hemoglobin in paroxysmal nocturnal hemoglobinuria. Transfusion 1997;37(3):291-7.
- 19** Mundee Y, Bigelow NC, et al. Simplified flow cytometric method for fetal hemoglobin containing red blood cells. Cytometry 2000;42(6):389-393.
- 20** Bromilow, IM,Duguid JK. Measurement of feto-maternal haemorrhage: a comparative study of three Kleihauer techniques and two flow cytometry methods. Clin & Lab Haematol, 1997;19(2):137-42.
- 21** Chen J, Bigelow N, et al. Proposed flow cytometric reference method for the determination of erythroid F-cell counts. Cytometry 2000;42(4):239-46.
- 22** Davis, BH, Davis, KT. Laboratory assessment of fetomaternal hemorrhage is improved using flow cytometry. Lab Med 2007;38:365-73.
- 23** Garner C, Tatoo T, Reittie JE et al. Genetic influences on F cells and other hematologic variables: a twin heritability study. Blood 2000;95: 342-46.
- 24** Maier-Redelsperger M, de Montalembert M, et al., Fetal hemoglobin and F-cell responses to long-term hydroxyurea treatment in young sickle cell patients. The French Study Group on Sickle Cell Disease. Blood 1998;91(12):4472-9.
- 25** Thein SL, Craig JE. Genetics of Hb F/F cell variance in adults and heterocellular hereditary persistence of fetal hemoglobin. Hemoglobin 1998;22(5-6):401-14
- 26** EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.
- 27** Corcoran D, Murphy D, Donnelly JC, Ainle FN. The prevalence of maternal F cells in a pregnant population and potential overestimation of foeto-maternal haemorrhage as a consequence. Blood Transfus. 2014 Oct;12(4):570-4
- 28** Pastoret C, Le Priol J, Fest T, Roussel Evaluation of FMH QuikQuant for the detection and quantification of fetomaternal hemorrhage. Cytometry B Clin Cytom. 2013 Jan- Feb;84(1):37-43

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Referencias

- 1** CLSI, Fetal red cell detection; approved guideline. CLSI (formerly NCCLS) Document H52-A, 2001.
- 2** Davis, BH. Diagnostic advances in defining erythropoietic abnormalities and red cell diseases. Seminars in Hematology, 2001;38:148-59.
- 3** Davis, B.H. Diagnostic utility of red cell flow cytometric analysis. Clin Lab Med 2001;21(4):829-40.
- 4** Sebring ES, Polesky HF. Fetomaternal haemorrhage: incidence, risk factors, time of occurrence, and clinical effects. Transfusion 1990;30:344-57.
- 5** Giacoia G.P. Severe fetomaternal hemorrhage: a review. Obstet & Gynecol Surv, 1997;52(6): p. 372-80.
- 6** Polesky, HF, Sebring ES, Evaluation of methods for detection and quantitation of fetal cells and their effect on RhIgG. American Journal of Clin Path, 1981;76:525-29.
- 7** Hartwell, EA. Use of Rh immune globulin: ASCP practice parameter. American Journal of Clin Path, 1998;110:281-92.
- 8** Lee D, Contreras M, Robson SC, et al. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. Transfusion Med 1999;9:93-7.
- 9** Kleihauer E, Braun H, Betke K. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. Klin Wochenschr 1957;35:637-38.
- 10** Duckett JR, Constantine G. The Kleihauer technique: an accurate method of quantifying fetomaternal haemorrhage? British Journal of Obstet & Gynaecol 1997;104: 845-6.
- 11** Emery CL, Morway LF , et al. The Kleihauer-Betke test. Clinical utility, indication, and correlation in patients with placental abruption and cocaine use. Arch Pathol Lab Med 1995;119(11):1032-7.
- 12** Davis BH, Olsen S, et al. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. Transfusion 1998;38(8):749-56.
- 13** Chen JC, Davis BH, et al. Multicenter clinical experience with flow cytometric method for fetomaternal hemorrhage detection. Cytometry 2002;50(6):285-90.
- 14** Lloyd-Evans P, Kumpel BM et al. Use of a directly conjugated monoclonal anti-D (BRAD-3) for quantification of fetomaternal hemorrhage by flow cytometry. Transfusion 1996;36(5):432-7.
- 15** Nance SJ, Nelson JM, et al. Quantitation of fetal-maternal hemorrhage by flow cytometry. A simple and accurate method. Am J Clin Pathol, 1989;91(3): 288-92.
- 16** Navenot JM, Merghoub T et al. New method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle-cell disease. Cytometry 1998;32(3):186-90.
- 17** Nelson M, Zarkos K et al. A flow-cytometric equivalent of the Kleihauer test. Vox Sanguinis 1998;75:234-41.
- 18** Navenot JM, Muller JY, et al. Expression of blood group i antigen and fetal hemoglobin in paroxysmal nocturnal hemoglobinuria. Transfusion 1997;37(3):291-7.
- 19** Mundee Y, Bigelow NC, et al. Simplified flow cytometric method for fetal hemoglobin containing red blood cells. Cytometry 2000;42(6):389-393.
- 20** Bromilow, IM,Duguid JK. Measurement of feto-maternal haemorrhage: a comparative study of three Kleihauer techniques and two flow cytometry methods. Clin & Lab Haematol, 1997;19(2):137-42.
- 21** Chen J, Bigelow N, et al. Proposed flow cytometric reference method for the determination of erythroid F-cell counts. Cytometry 2000;42(4):239-46.
- 22** Davis, BH, Davis, KT. Laboratory assessment of fetomaternal hemorrhage is improved using flow cytometry. Lab Med 2007;38:365-73.
- 23** Garner C, Tatu T, Reittie JE et al. Genetic influences on F cells and other hematologic variables: a twin heritability study. Blood 2000;95: 342-46.
- 24** Maier-Redelsperger M, de Montalembert M, et al., Fetal hemoglobin and F-cell responses to long-term hydroxyurea treatment in young sickle cell patients. The French Study Group on Sickle Cell Disease. Blood 1998;91(12):4472-9.
- 25** Thein SL, Craig JE. Genetics of Hb F/F cell variance in adults and heterocellular hereditary persistence of fetal hemoglobin. Hemoglobin 1998;22(5-6):401-14
- 26** EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.
- 27** Corcoran D, Murphy D, Donnelly JC, Ainle FN. The prevalence of maternal F cells in a pregnant population and potential overestimation of foeto-maternal haemorrhage as a consequence. Blood Transfus. 2014 Oct;12(4):570-4
- 28** Pastoret C, Le Priol J, Fest T, Roussel Evaluation of FMH QuikQuant for the detection and quantification of fetomaternal hemorrhage. Cytometry B Clin Cytom. 2013 Jan- Feb;84(1):37-43

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References

- 1 CLSI, Fetal red cell detection; approved guideline. CLSI (formerly NCCLS) Document H52-A, 2001.
- 2 Davis, BH. Diagnostic advances in defining erythropoietic abnormalities and red cell diseases. Seminars in Hematology, 2001;38:148-59.
- 3 Davis, B.H. Diagnostic utility of red cell flow cytometric analysis. Clin Lab Med 2001;21(4):829-40.
- 4 Sebring ES, Polesky HF. Fetomaternal haemorrhage: incidence, risk factors, time of occurrence, and clinical effects. Transfusion 1990;30:344-57.
- 5 Giacoia G.P. Severe fetomaternal hemorrhage: a review. Obstet & Gynecol Surv, 1997;52(6): p. 372-80.
- 6 Polesky, HF, Sebring ES, Evaluation of methods for detection and quantitation of fetal cells and their effect on RhIgG. American Journal of Clin Path, 1981;76:525-29.
- 7 Hartwell, EA. Use of Rh immune globulin: ASCP practice parameter. American Journal of Clin Path, 1998;110:281-92.
- 8 Lee D, Contreras M, Robson SC, et al. Recommendations for the use of anti-D immunoglobulin for Rh prophylaxis. Transfusion Med 1999;9:93-7.
- 9 Kleihauer E, Braun H, Betke K. Demonstration of fetal hemoglobin in erythrocytes of a blood smear. Klin Wochenschr 1957;35:637-38.
- 10 Duckett JR, Constantine G. The Kleihauer technique: an accurate method of quantifying fetomaternal haemorrhage? British Journal of Obstet & Gynaecol 1997;104: 845-6.
- 11 Emery CL, Morway LF , et al.The Kleihauer-Betke test. Clinical utility, indication, and correlation in patients with placental abruption and cocaine use. Arch Pathol Lab Med 1995;119(11):1032-7.
- 12 Davis BH, Olsen S, et al. Detection of fetal red cells in fetomaternal hemorrhage using a fetal hemoglobin monoclonal antibody by flow cytometry. Transfusion 1998;38(8):749-56.
- 13 Chen JC, Davis BH, et al. Multicenter clinical experience with flow cytometric method for fetomaternal hemorrhage detection. Cytometry 2002;50(6):285-90.
- 14 Lloyd-Evans P, Kumpel BM et al. Use of a directly conjugated monoclonal anti-D (BRAD-3) for quantification of fetomaternal hemorrhage by flow cytometry. Transfusion 1996;36(5):432-7.
- 15 Nance SJ, Nelson JM, et al. Quantitation of fetal-maternal hemorrhage by flow cytometry. A simple and accurate method. Am J Clin Pathol, 1989;91(3): 288-92.
- 16 Navenot JM, Merghoub T et al. New method for quantitative determination of fetal hemoglobin-containing red blood cells by flow cytometry: application to sickle-cell disease. Cytometry 1998;32(3):186-90.
- 17 Nelson M, Zarkos K et al. A flow-cytometric equivalent of the Kleihauer test. Vox Sanguinis 1998;75:234-41.
- 18 Navenot JM, Muller JY, et al. Expression of blood group i antigen and fetal hemoglobin in paroxysmal nocturnal hemoglobinuria. Transfusion 1997;37(3):291-7.
- 19 Mundee Y, Bigelow NC, et al. Simplified flow cytometric method for fetal hemoglobin containing red blood cells. Cytometry 2000;42(6):389-393.
- 20 Bromilow, IM,Duguid JK. Measurement of feto-maternal haemorrhage: a comparative study of three Kleihauer techniques and two flow cytometry methods. Clin & Lab Haematol, 1997;19(2):137-42.
- 21 Chen J, Bigelow N, et al. Proposed flow cytometric reference method for the determination of erythroid F-cell counts. Cytometry 2000;42(4):239-46.
- 22 Davis, BH, Davis, KT. Laboratory assessment of fetomaternal hemorrhage is improved using flow cytometry. Lab Med 2007;38:365-73.
- 23 Garner C, Tatu T, Reittie JE et al. Genetic influences on F cells and other hematologic variables: a twin heritability study. Blood 2000;95: 342-46.
- 24 Maier-Redelsperger M, de Montalembert M, et al., Fetal hemoglobin and F-cell responses to long-term hydroxyurea treatment in young sickle cell patients. The French Study Group on Sickle Cell Disease. Blood 1998;91(12):4472-9.
- 25 Thein SL, Craig JE. Genetics of Hb F/F cell variance in adults and heterocellular hereditary persistence of fetal hemoglobin. Hemoglobin 1998;22(5-6):401-14
- 26 EN ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.
- 27 Corcoran D, Murphy D, Donnelly JC, Ainle FN. The prevalence of maternal F cells in a pregnant population and potential overestimation of foeto-maternal haemorrhage as a consequence. Blood Transfus. 2014 Oct;12(4):570-4
- 28 Pastoret C, Le Priol J, Fest T, Roussel Evaluation of FMH QuikQuant for the detection and quantification of fetomaternal hemorrhage. Cytometry B Clin Cytom. 2013 Jan- Feb;84(1):37-43

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Tilauskia: orders@iqproducts.nl

www.iqproducts.nl

©2023 - IQ Products bv. All rights reserved. No parts of these works may be reproduced in any form without permission in writing.

Current version + release date	Version 3 29-03-2023 (DD-MM-YYYY)
Previous version	Version 2
Changes	<ul style="list-style-type: none"> • Removed sentence about IVD registration on page 1, 10, 18 and 26. • Added Switzerland and Peru as IVD registered countries
Justification	<ul style="list-style-type: none"> • This sentence should have been deleted with expanding the part about regulatory status in the previous version. • Complete information.
Previous version + release date	Version 2 20-02-2023 (DD-MM-YYYY)
Previous version	Version 1
Changes	Added Australia and United Kingdom as countries where the FMH QuikQuant has been registered as IVD product.
Justification	This information was missing and is useful for users in these countries.
Previous version + release date	Version 1 01-09-2021 (DD-MM-YYYY)
Previous version	N/A
Changes	1. Created this online version with translations of the international booklet which is printed and included with every kit. The online version includes Portuguese, Spanish and Finnish, whereas the international booklet includes English, French and German.
Justification	1. Environmental reasons.