CMV Brite™ Turbo Kit

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1. Van der Bij, W., Schirm, J., Torensma, R., Van Son, W.J., Tegzess, A.M., The, T.H. (1988) Comparison between viremia and antigenemia for detection of cytomegalovirus in blood. J. Clin. Microbiology **26**, 2531 - 2535.
2. Grefte, J.M.M., Van der Gun, B.T.F., Schmolke, S., Van der Giessen, M., Van Son, W.J. and The, T.H. (1992). The lower matrix protein pp65 is the principal viral antigen present in peripheral blood leukocytes during an active cytomegalovirus infection. J.Gen.Virol. **73**, 2923-2932.
3. Gerna, G., Revello, M.G., Percivalle, E., Zavattoni, M., Parea, M., Battaglia, M. (1990). Quantification of human cytomegalovirus viremia by using monoclonal antibodies to different viral proteins. J. Clin. Microbiology **28**, 2681 - 2688.
4. Revello, M.G., Percivalle, E., Di Matteo, A., Morini, F., Gerna, G. (1992). Nuclear expression of the lower matrix protein of human cytomegalovirus in peripheral blood leukocytes of immunocompromised viremic patients. J. Gen Virol **73**, 437 - 442.
5. Gerna, G., Revello, M.G., Percivalle, E., Morini, F. (1992). Comparison of different immunostaining techniques and monoclonal antibodies to the lower matrix phosphoprotein (pp65) for optimal quantitation of human cytomegalovirus antigenemia. J. Clin. Microbiology **30**, 1232 - 1237.
6. Ho, S.K.N., Lo, C-Y., Cheng, I.K.P. and Chan, T-M., (1998) Rapid cytomegalovirus pp65 antigenemia assay by direct Erythrocyte Lysis and Immunofluorescence Staining. J. Clin.Microbiology. **36**, 638-640.
7. Boeckh, M., Woogerd, P.M., Stevens-Ayers, T., Ray, C.G., and Bowden, R.A. (1994). Factors influencing detection of quantitative Cytomegalovirus antigenemia. J.Clin.Microbiology **32**, 832-834.
8. Landry, M.L., Ferguson D., Cohen, S., Huber, K., and Wetherill, P. (1995). Effect of delayed specimen processing on Cytomegalovirus antigenemia test results. J.Clin.Microbiology **33**, 257-259.
9. Landry, M.L., Ferguson, D., Stevens-Ayers, T., De Jonge, M.W.A., and Boeckh, M. (1996). Evaluation of CMV Brite™ Kit for detection of Cytomegalovirus pp65 antigenemia in peripheral blood leukocytes by immunofluorescence. J.Clin.Microbiology **34**, 1337-1339.
10. ISO 15223-1:2021 Ιατρικές συσκευές – Σύμβολα που πρέπει να χρησιμοποιούνται με τις πληροφορίες που πρέπει να παρέχονται από τον

κατασκευαστή - Μέρος 1: Γενικές απαιτήσεις.Niubò, J., Pérez, J.L., Carvajal, A., Ardanuy, C. and Martín, R. (1994). Effect of delayed processing of blood samples on performance of cytomegalovirus antigenemia assay. J.Clin.Microbiology **32**, 1119-1120.

11. ISO 15223-1:2021/Amd 1:2025 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements - Amendment 1: Addition of defined term for authorized representative and modified EC REP symbol to not be country or region specific.

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- Van der Bij, W., Schirm, J., Torensma, R., Van Son, W.J., Tegzess, A.M., The, T.H. (1988) Comparison between viremia and antigenemia for detection of cytomegalovirus in blood. J. Clin. Microbiology 26, 2531 - 2535.
- Grefte, J.M.M., Van der Gun, B.T.F., Schmolke, S., Van der Giessen, M., Van Son, W.J. and The, T.H. (1992). The lower matrix protein pp65 is the principal viral antigen present in peripheral blood leukocytes during an active cytomegalovirus infection. J.Gen.Virol. 73, 2923-2932.
- Gerna, G., Revello, M.G., Percivalle, E., Zavattoni, M., Parea, M., Battaglia, M. (1990). Quantification of human cytomegalovirus viremia by using monoclonal antibodies to different viral proteins. J. Clin. Microbiology 28, 2681 - 2688.
- Revello, M.G., Percivalle, E., Di Matteo, A., Morini, F., Gerna, G. (1992). Nuclear expression of the lower matrix protein of human cytomegalovirus in peripheral blood leukocytes of immunocompromised viremic patients. J. Gen Virol 73, 437 - 442.
- Gerna, G., Revello, M.G., Percivalle, E., Morini, F. (1992). Comparison of different immunostaining techniques and monoclonal antibodies to the lower matrix phosphoprotein (pp65) for optimal quantitation of human cytomegalovirus antigenemia. J. Clin. Microbiology 30, 1232 - 1237.
- Ho, S.K.N., Lo, C-Y., Cheng, I.K.P. and Chan, T-M., (1998) Rapid cytomegalovirus pp65 antigenemia assay by direct Erythrocyte Lysis and Immunofluorescence Staining. J. Clin.Microbiology. 36, 638-640.
- Boeckh, M., Woogerd, P.M., Stevens-Ayers, T., Ray, C.G., and Bowden, R.A. (1994). Factors influencing detection of quantitative Cytomegalovirus antigenemia. J.Clin.Microbiology 32, 832-834.
- Landry, M.L., Ferguson, D., Cohen, S., Huber, K., and Wetherill, P. (1995). Effect of delayed specimen processing on Cytomegalovirus antigenemia test results. J.Clin.Microbiology 33, 257-259.
- Landry, M.L., Ferguson, D., Stevens-Ayers, T., De Jonge, M.W.A., and Boeckh, M. (1996). Evaluation of CMV Brite™ Kit for detection of Cytomegalovirus pp65 antigenemia in peripheral blood leukocytes by immunofluorescence. J.Clin.Microbiology 34, 1337-1339.
- NEN EN ISO 15223-1 Dispositivos médicos - Símbolos que se utilizarán con etiquetas de dispositivos médicos, el etiquetado y la información a suministrar - Parte 1: Requisitos generales.
- Niubò, J., Pérez, J.L., Carvajal, A., Ardanuy, C. and Martín, R. (1994). Effect of delayed processing of blood samples on performance of cytomegalovirus antigenemia assay. J.Clin.Microbiology 32, 1119-1120.
- ISO 15223-1:2021/Amd 1:2025 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements - Amendment 1: Addition of defined term for authorized representative and modified EC REP symbol to not be country or region specific.

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Bibliografia

1. Van der Bij, W., Schirm, J., Torensma, R., Van Son, W.J., Tegzess, A.M., The, T.H. (1988) Comparison between viremia and antigenemia for detection of cytomegalovirus in blood. J. Clin. Microbiology 26, 2531 - 2535.
2. Grefte, J.M.M., Van der Gun, B.T.F., Schmolke, S., Van der Giessen, M., Van Son, W.J. and The, T.H. (1992). The lower matrix protein pp65 is the principal viral antigen present in peripheral blood leukocytes during an active cytomegalovirus infection. J.Gen.Virol. 73, 2923-2932.
3. Gerna, G., Revello, M.G., Percivalle, E., Zavattoni, M., Parea, M., Battaglia, M. (1990). Quantification of human cytomegalovirus viremia by using monoclonal antibodies to different viral proteins. J. Clin. Microbiology 28, 2681 - 2688.
4. Revello, M.G., Percivalle, E., Di Matteo, A., Morini, F., Gerna, G. (1992). Nuclear expression of the lower matrix protein of human cytomegalovirus in peripheral blood leukocytes of immunocompromised viremic patients. J. Gen Virol 73, 437 - 442.
5. Gerna, G., Revello, M.G., Percivalle, E., Morini, F. (1992). Comparison of different immunostaining techniques and monoclonal antibodies to the lower matrix phosphoprotein (pp65) for optimal quantitation of human cytomegalovirus antigenemia. J. Clin. Microbiology 30, 1232 - 1237.
6. Ho, S.K.N., Lo, C-Y., Cheng, I.K.P. and Chan, T-M., (1998) Rapid cytomegalovirus pp65 antigenemia assay by direct Erythrocyte Lysis and Immunofluorescence Staining. J. Clin. Microbiology. 36, 638-640.
7. Boeckh, M., Woogerd, P.M., Stevens-Ayers, T., Ray, C.G., and Bowden, R.A. (1994). Factors influencing detection of quantitative Cytomegalovirus antigenemia. J.Clin.Microbiology 32, 832-834.
8. Landry, M.L., Ferguson, D., Cohen, S., Huber, K., and Wetherill, P. (1995). Effect of delayed specimen processing on Cytomegalovirus antigenemia test results. J.Clin.Microbiology 33, 257-259.
9. Landry, M.L., Ferguson, D., Stevens-Ayers, T., De Jonge, M.W.A., and Boeckh, M. (1996). Evaluation of CMV Brite TM Kit for detection of Cytomegalovirus pp65 antigenemia in peripheral blood leukocytes by immunofluorescence. J.Clin.Microbiology 34, 1337-1339.
10. ISO 15223-1:2021 Medical devices – Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements.
11. Niubó, J., Pérez, J.L., Carvajal, A., Ardanuy, C. and Martín, R. (1994). Effect of delayed processing of blood samples on performance of cytomegalovirus antigenemia assay. J.Clin.Microbiology 32, 1119-1120.
12. ISO 15223-1:2021/Amd 1:2025 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements - Amendment 1: Addition of defined term for authorized representative and modified EC REP symbol to not be country or region specific.

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Current version + release date	Version 6 21-01-2026 (DD-MM-YYYY)
Previous version	Version 5
Changes	1. added ISO 15223-1:2021/Amd 1:2025 to bibliography 2. replaced symbol EC-REP for XX-REP
Justification	1. and 2. Publication of ISO 15223-1:2021/Amd 1:2025
Current version + release date	Version 5 01-11-2025 (DD-MM-YYYY)
Previous version	Version 4
Changes	1. The CE symbol is no longer shown in italics 2. hazard symbols of reagent B have been updated 3. re-arranged information on front page 4. Lay-out of table 'Kit content' has been updated 5. the word 'non-automated' was added to the intended use. 6. removed CE from table with Explanation of used symbols
Justification	1. In response to a non-conformity raised by the notified body. 2. In response to the CLP analysis performed for the IVDR Technical Documentation. 3. Increase readability 4. Remove non-existing symbols 5. IVDR requirements for components of an intended use. 6. CE is not a symbol, but an abbreviation
Current version + release date	Version 4 February 2025
Previous version	Version 3
Changes	1. Sentence "This product is not FDA cleared (approved) for use in testing (i.e. screening) of blood or plasma donors." Was moved from sections Intended use to section Regulatory status. Applicability of the sentence has been broadened by removing the FDA part. 2. Added Australia, Peru, Saudi-Arabia and Singapore as country in which this product is registered as IVD.
Justification	1. Regulatory information is not supposed to be mentioned in the Intended Use. This product is not intended to be used for screening in the entire world, meaning that the specification for the FDA is redundant. 2. Registration in these countries is complete.
Current version + release date	Version 3 December 2024
Previous version	Version 2
Changes	1. Moved heading 'Regulatory status'. 2. Correct CI of CMV Brite and added sensitivity and specificity (including CI) for CMV Brite Turbo.
Justification	1. Making this information more visible. 2. Described CI of CMV Brite could not be reproduced. Information on CMV Brite Turbo missed, but was requested for registration of this product in Singapore.

Current version + release date	Version 2 March 2024
Previous version	Version 1
Changes	GTIN and UDI-DI code on front page have been corrected.
Justification	Wrong codes were included.

Current version + release date	Version 1 April 2023
Previous version	N/A
Changes	1. Created this online version with translations of the international booklet which is printed and included with every kit. The online version includes Greek, Spanish and Italian, whereas the international booklet includes English and Portuguese.
Justification	1. Environmental reasons.

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©2026 - IQ Products BV. Όλα τα δικαιώματα διατηρούνται. Κανένα από τα τμήματα των έργων αυτών δεν μπορεί να αναπαραχθεί σε οποιαδήποτε μορφή χωρίς δικαιώματα.