

PRODUCT INFORMATION SHEET

Mouse Mo	e Monoclonal Antibody to HAUSP							
	Format: PURE	RUO	REF	IQP-1440P	Volume:	0.1ml		
	Description							
Clone	5F11							
Isotype	Mouse IgG1							
Specificity	HAUSP							
Alternative names	TEF1; HAUSP; USP7							
Species	Human Purified recombinant fragment of human HAUSP expressed in E. Coli.							
Immunogen		-						
Mw	128kDa							
Format	Ascitic fluid containing 0.03% sodium azide.							
Summary	USP7 or HAUSP is a ubiquitin specific protease or a deubiquitylating enzyme that cleaves ubiquitin from its substrates. Since ubiquitylation (polyubiquitination) is most commonly associated with the stability and degradation of cellular proteins, HAUSP acitivity generally stabilizes its substrate proteins. HAUSP is most popularly known as a direct antagonist of Mdm2, the E3 ubiquitin ligase for the tumor suppressor protein, p53.Normally, p53 levels are kept low in part due to Mdm2- mediated ubiquitylation and degradation of p53. Interestingly, in response to oncogenic insults, HAUSP can deubiquitinate p53 and protect p53 from Mdm2-mediated degradation, indicating that it may possess a tumor suppressor function for the immediate stabilization of p53 in response to stress. Another important role of HAUSP function involves the oncogenic stabilization of p53. Oncogenes such as Myc and E1A are thought to activate p53 through a p19 alternative reading frame (p19ARF, also called ARF)-dependent pathway, although some evidence suggests ARF is no essential in this process. An intriguing possibility is that HAUSP provides an alternative pathway for safeguarding the cell against oncogenic insults.							
Applications	Western Bloting: 1/500 - 1/2000.ELISA: Propose dilution 1/10000.Not yet tested in other applications.Determining optimal working dilutions by titration test.							
Limitations	 Conjugates with brighter fluorochromes, like PE and APC, will have a greater separation tha those with dyes like FITC and CyQ. When populations overlap, the percentage of positive ceusing a selected marker can be affected by the choice of fluorescent label. Use of monoclonal antibodies in patient treatment can interfere with antigen target recognit by this reagent. This should be taken into account when samples are analyzed from patients treated in this fashion. IQ Products has not characterized the effect of the presence of therapeutic antibodies on the performance of this reagent. Reagents can be used in different combinations, therefore laboratories need to become familiar performance characteristics of each antibody in relation with the combined markers normal and abnormal samples. Reagent data performance is based on EDTA-treated blood. Reagent performance can be affected by the use of other anticoagulants. 							

⚠ 💩 🔏 巻 🎴

Handling and Storage

Antibodies are supplied in 0.01 M sodium phosphate, 0.15 M NaCl; pH 7.3, 0.03% sodiumazide (NaN3) or as ascetic fluid containing 0.03% sodiumazide. Store the vials at 2-8 °C for a maximum of 2 weeks and store at -20°C for longer term storage. Monoclonal antibodies should be protected from prolonged exposure to light when conjugated with fluorochromes. Reagents are stable for the period shown on the vial label when stored properly.

Warranty Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. IQ Products is not liable for property damage, personal injury, or economic loss caused by the product.

Characterization

- To ensure consistently high-quality reagents, each batch of monoclonal antibody is tested for conformance with characteristics of a standard reagent. Representative flow cytometric data is included in this data sheet.
- **Warning** All products contain sodiumazide. This chemical is poisonous and hazardous. Handling should be done by trained staff only.

References 1. Cell Death Differ. 2007 Jul;14(7):1350-60.

- 2. Cancer Cell. 2007 Oct;12(4):342-54.
- 3. Blood. 2009 Apr 2;113(14):3264-75.

Explanation of used symbols

	Consult instructions for use
REF V	Catalogue number
V	Sufficient for
Å	
	Caution, consult accompanying document
*	Keep away from (sun)light
æ	Biological risks
	Temperature limitation (°C)
RUO	For Research Use Only
RUO LOT	Batch code
2	Use by yyyy-mm-dd
	Manufacturer
EC REP	Authorized Representative in the European Community
CE	Conformité Européenne (European Conformity)

		Label - tandem	Ex -max (nm)	Em -max (nm)
Р	PURE	purified material	-	-
F	FITC	FITC	488	519
R	R-PE	PE	488, 532	578
С	CyQ	PE-Cy5.18	488, 532	667
Α	APC		595, 633, 635, 647	660
PC	PerCP		488, 532	678
PCC	PerCP-Cy5.5		488, 532	695

IQ Products BV

Rozenburglaan 13a 9727 DL Groningen, The Netherlands

+31 (0)50 57 57 000

- -+31 (0)50 57 57 002
- Technical <u>marketing@iqproducts.nl</u>
- Orders <u>orders@iqproducts.nl</u>
- <u>www.iqproducts.nl</u>