

PRODUCT INFORMATION SHEET

Mouse Monoclonal Antibody to HK2

Format: PURE RUO REF IQP-1451P Volume: 0.1ml

Description

Clone 3D3

Isotype Mouse IgG1

Specificity HK2

Alternative names

HKII; HXK2; DKFZp686M1669; HK2

Species Human

Immunogen Purified recombinant fragment of human HK2 expressed in E. Coli.

Mw 102kDa

Format Ascitic fluid containing 0.03% sodium azide.

Summary

The hexokinases utilize Mg-ATP as a phosphoryl donor to catalyze the first step of intracellular glucose metabolism, the conversion of glucose to glucose- 6-phosphate. Four hexokinase isoenzymes have been identified, including hexokinase I (HXK I), hexokinase II (HXK II), hexokinase III (HXK III) and hexokinase IV (HXK IV, also designated glucokinase or GCK). Hexokinases I-III each contain an N-terminal cluster of hydrophobic amino acids. Glucokinase lacks the N-terminal hydrophobic cluster. The hydrophobic cluster is thought to be necessary for membrane binding. This is substantiated by the finding that glucokinase has lower affinity for glucose than do the other hexokinases .Hexokinase 2 is the predominant hexokinase isozyme expressed in insulin-responsive tissues such as skeletal muscle. Expression of this gene is insulin-responsive, and studies in rat suggest that it is involved in the increased rate of glycolysis seen in rapidly growing cancer cells.

Applications

Western Bloting: 1/500 - 1/2000.Immunohistochemistry: 1/200 - 1/1000.Flow cytometry: 1/200 - 1/400.ELISA: Propose dilution 1/10000.Not yet tested in other applications.Determining optimal working dilutions by titration test.

Limitations

- 1. Conjugates with brighter fluorochromes, like PE and APC, will have a greater separation than those with dyes like FITC and CyQ. When populations overlap, the percentage of positive cells using a selected marker can be affected by the choice of fluorescent label.
- Use of monoclonal antibodies in patient treatment can interfere with antigen target recognition by this reagent. This should be taken into account when samples are analyzed from patients treated in this fashion. IQ Products has not characterized the effect of the presence of therapeutic antibodies on the performance of this reagent.
- 3. Reagents can be used in different combinations, therefore laboratories need to become familiar performance characteristics of each antibody in relation with the combined markers in normal and abnormal samples.
- 4. Reagent data performance is based on EDTA-treated blood. Reagent performance can be affected by the use of other anticoagulants.

IOP-1451 - HK2 Version 1

△ & / * □

Handling and Storage

Antibodies are supplied in 0.01 M sodium phosphate, 0.15 M NaCl; pH 7.3, 0.03% sodiumazide (NaN3) or as ascetic fluid containing 0.03% sodiumazide. Store the vials at 2-8 °C for a maximum of 2 weeks and store at -20°C for longer term storage. Monoclonal antibodies should be protected from prolonged exposure to light when conjugated with fluorochromes. Reagents are stable for the period shown on the vial label when stored properly.

Warranty

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. IQ Products is not liable for property damage, personal injury, or economic loss caused by the product.

Characterization

To ensure consistently high-quality reagents, each batch of monoclonal antibody is tested for conformance with characteristics of a standard reagent. Representative flow cytometric data is included in this data sheet.

Warning

All products contain sodiumazide. This chemical is poisonous and hazardous. Handling should be done by trained staff only.

References 1. Cell. 2006 May 19;125(4):801-14.

Cancer Sci. 2008 Feb;99(2):260-6.
Med Oncol. 2009;26(3):303-8.

Explanation of used symbols

Consult instructions for use REF Catalogue number Sufficient for \triangle Caution, consult accompanying document 巻 Keep away from (sun)light 8 Biological risks Temperature limitation (°C) right fluorescence RUO For Research Use Only LOT Batch code Use by yyyy-mm-dd Manufacturer **ECREP** Authorized Representative in the European Community

P F R C	PURE FITC R-PE CyQ APC	Label - tandem purified material FITC PE PE-Cy5.18	Ex -max (nm) - 488 488, 532 488, 532 595, 633, 635, 647	Em -max (nm) - 519 578 667 660
PC	PerCP		488, 532	678
PCC	PerCP-Cy5.5		488, 532	695



CE

IQ Products BV

Rozenburglaan 13a 9727 DL Groningen, The Netherlands

** +31 (0)50 57 57 000 ** +31 (0)50 57 57 002

↑ Technical marketing@iqproducts.nl
↑ Orders orders@iqproducts.nl

Conformité Européenne (European Conformity)

<u>www.iqproducts.nl</u>

IOP-1451 - HK2 Version 1