

PRODUCT INFORMATION SHEET

Monoclonal antibodies detecting human antigens

Anti-HLA-G

PURE RUO REF IQP-634P ▼ 100 test

RUO For Research Use Only

Description

Clone MEM-G/2
Isotype Murine IgG1

The antibody MEM-G/2 recognizes the free heavy chain of all the HLA-G isoforms. HLA-G belongs to the MHC Class I molecules (MHC Class Ib; nonclassical) and it is expressed on the surface of

trophoblast cells.

Species Human

Immunogen Bacterially expressed recombinant HLA-G heavy chain (denatured).

Applications IHC(P), WB, IHC(F). Determining optimal working dilutions by titration test.

Limitations

1. Conjugates with brighter fluorochromes, like PE and APC, will have a greater separation than those with dyes like FITC and CyQ. When populations overlap, the percentage of positive cells using a selected marker can be affected by the choice of fluorescent label.

2. Use of monoclonal antibodies in patient treatment can interfere with antigen target recognition by this reagent. This should be taken into account when samples are analyzed from patients treated in this fashion. IQ Products has not characterized the effect of the presence of therapeutic antibodies on the performance of this reagent.

3. Reagents can be used in different combinations, therefore laboratories need to become familiar performance characteristics of each antibody in relation with the combined markers in normal and abnormal samples.

4. Reagent performance can be affected by the use of anticoagulants.

△ ♦ / * 2

Handling and Storage

Antibodies are supplied in phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4. Store the vials at 2-8°C. Monoclonal antibodies should be protected from prolonged exposure to light when conjugated with fluorochromes. Reagents are stable for the period shown on the vial label when stored properly.

Warranty

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. IQ Products is not liable for property damage, personal injury, or economic loss caused by the product.

Characterization

To ensure consistently high-quality reagents, each batch of monoclonal antibody is tested for conformance with characteristics of a standard reagent.

Warning

All products contain sodiumazide. This chemical is poisonous and hazardous. Handling should be done by trained staff only.

References

 Menier C, Saez B, Horejsi V, Martinozzi S, Krawice-Radanne I, Bruel S, LeDanff C, Reboul M, Hilgert I, Rabreau M, Larrad ML, Pla M, Carosella ED, Rouas-Freiss N: Characterization of monoclonal antibodies recognizing HLA-G or HLA-E: new tools to analyze the expression of nonclassical HLA class I molecules. Hum Immunol. 2003 64(3):315-26.

- 2. Creput C, Durrbach A, Menier C, Guettier C, Samuel D, Dausset J, Charpentier B, Carosella ED, Rouas-Freiss N. Human leukocyte antigen-G (HLA-G) expression in biliary epithelial cells is associated with allograft acceptance in liver-kidney transplantation. J Hepatol. 2003 39(4):587-94.
- 3. Boyson JE, Erskine R, Whitman MC, Chiu M, Lau JM, Koopman LA, Valter MM, Angelisova P, Horejsi V, Strominger JL: Disulfide bond-mediated dimerization of HLA-G on the cell surface. Proc Natl Acad Sci U S A. 2002 Dec 10;99(25):16180-5.
- 4. Lazana I, Zoudiari A, Kokkinou D, Themeli M, Liga M, Papadaki H, Papachristou D, Spyridonidis A: Identification of a novel HLA-G+ regulatory population in blood: expansion after allogeneic transplantation and de novo HLA-G expression at graft-versus-host disease sites. Haematologica. 2012 Sep;97(9):1338-47
- 5. Poláková K, Bandzuchová E, Hofmeister V, Weiss EH, Hutter H, Russ G: Binding analysis of HLA-G specific antibodies to hematopoietic cells isolated from leukemia patients. Neoplasma. 2003;50(5):331-8.
- 6. Le Discorde M, Moreau P, Sabatier P, Legeais JM, Carosella ED: Expression of HLA-G in human cornea, an immune-privileged tissue. Hum Immunol. 2003 Nov;64(11):1039-44.

Explanation of used symbols

(li	Consult instructions for use
REF	Catalogue number
\$	Sufficient for
\triangle	Caution, consult accompanying document
*	Keep away from (sun)light
€	Biological risks
*	Temperature limitation (°C)
RUO	For Research Use Only
LOT	Batch code
\square	Use by yyyy-mm-dd
	Manufacturer

		Label - tandem	Ex -max (nm)	Em -max (nm)
P	PURE	purified material	-	-
F	FITC	FITC	488	519
R	R-PE	PE	488, 532	578
C	CyQ	PE-Cy5.18	488, 532	667
Α	APC		595, 633, 635, 647	660
PC	PerCP		488, 532	678
PCC	PerCP-Cy5.5		488, 532	695

nt fluorescence



IQ Products BV

Rozenburglaan 13a 9727 DL Groningen, The Netherlands

- ♣ +31 (0)50 57 57 002
- Technical <u>marketing@iqproducts.nl</u>
- Orders orders@igproducts.nl
- <u>www.iqproducts.nl</u>