

PRODUCT INFORMATION SHEET

Mouse Monoclonal Antibody to HDAC3

Format: PURE RUO REF IQP-1218P Volume: 0.1ml

Description Ti

Clone 3A7B5

Isotype Mouse IgG1 HDAC3 **Specificity**

Alternative names

HD3; RPD3; RPD3-2

Species Human

Immunogen Purified recombinant fragment of HDAC3 (aa224-428) expressed in E. Coli.

48.8kDa

Format Ascitic fluid containing 0.03% sodium azide.

Summary

HDAC3: histone deacetylase 3, also known as HD3, RPD3, RPD3-2. Entrez Protein NC_000005. Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to the histone deacetylase/acuc/apha family. It has histone deacetylase activity and represses transcription when tethered to a promoter. It may participate in the regulation of transcription through its binding with the zinc-finger transcription factor YY1. This protein can also down-regulate p53 function and thus modulate cell growth and apoptosis. This gene is regarded as a potential tumor suppressor gene.

Applications Western Bloting: 1/500 - 1/2000.Immunohistochemistry: 1/200 - 1/1000.ELISA: Propose dilution 1/10000. Not yet tested in other applications. Determining optimal working dilutions by titration test.

Limitations

- 1. Conjugates with brighter fluorochromes, like PE and APC, will have a greater separation than those with dyes like FITC and CyQ. When populations overlap, the percentage of positive cells using a selected marker can be affected by the choice of fluorescent label.
- 2. Use of monoclonal antibodies in patient treatment can interfere with antigen target recognition by this reagent. This should be taken into account when samples are analyzed from patients treated in this fashion. IQ Products has not characterized the effect of the presence of therapeutic antibodies on the performance of this reagent.
- 3. Reagents can be used in different combinations, therefore laboratories need to become familiar performance characteristics of each antibody in relation with the combined markers in normal and abnormal samples.
- 4. Reagent data performance is based on EDTA-treated blood. Reagent performance can be affected by the use of other anticoagulants.

IQP-1218 - HDAC3 Version 1

▲ ♦ / * 2

Handling and Storage

Antibodies are supplied in 0.01 M sodium phosphate, 0.15 M NaCl; pH 7.3, 0.03% sodiumazide (NaN3) or as ascetic fluid containing 0.03% sodiumazide. Store the vials at 2-8 °C for a maximum of 2 weeks and store at -20°C for longer term storage. Monoclonal antibodies should be protected from prolonged exposure to light when conjugated with fluorochromes. Reagents are stable for the period shown on the vial label when stored properly.

Warranty

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. IQ Products is not liable for property damage, personal injury, or economic loss caused by the product.

Characterization

To ensure consistently high-quality reagents, each batch of monoclonal antibody is tested for conformance with characteristics of a standard reagent. Representative flow cytometric data is included in this data sheet.

Warning

巻

8

LOT

All products contain sodiumazide. This chemical is poisonous and hazardous. Handling should be done by trained staff only.

nt fluorescence

References 1. Folia Biol (Praha). 2006;52(1-2):21-33.

- 2. Mol Cell Biol. 2007 Feb;27(4):1280-95.
- 3. Mol Carcinog. 2008 Feb;47(2):137-47.

Explanation of used symbols



Sufficient for

 \triangle Caution, consult accompanying document

Keep away from (sun)light

Biological risks

Temperature limitation (°C) RUO For Research Use Only

Batch code

Use by yyyy-mm-dd

Manufacturer **ECREP**

Authorized Representative in the European Community Conformité Européenne (European Conformity)

		Label - tandem	Ex -max (nm)	Em -max (nm)
Р	PURE	purified material	-	-
F	FITC	FITC	488	519
R	R-PE	PE	488, 532	578
С	CyQ	PE-Cy5.18	488, 532	667
Α	APC	•	595, 633, 635, 647	660
PC	PerCP		488, 532	678
PCC	PerCP-Cy5.5		488, 532	695



IQ Products BV

Rozenburglaan 13a 9727 DL Groningen, The Netherlands

+31 (0)50 57 57 000 +31 (0)50 57 57 002

Technical marketing@igproducts.nl Orders orders@igproducts.nl

www.igproducts.nl

IQP-1218 - HDAC3 Version 1